

Sulfur (S) 1000 mg/L in H₂O, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003007
Unit 100 mL
Deliverydetails No Dangerous Good

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--|
| | Sulfur (S) | [7704-34-9] | 1000 | mg/L | | (NH ₄) ₂ SO ₄ , 99.999% |