

Tin (Sn) 1000 mg/L in 1% HNO3 0.5% HF, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003011
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) III E2

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	1000	mg/L		Sn, 99.999 9%