

Titanium (Ti) 1000 mg/L in diluted HNO₃/HF tr, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003014
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Titanium (Ti)	[7440-32-6]	1000	mg/L		(NH ₄) ₂ TiF ₆ ,, 99.998%