

**Titanium (Ti) 1000 mg/L in diluted HNO3/HF tr, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4003014  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source   |
|------------------|-------------------|-------------|---------------------|------|--------|--|
|                  | Titanium (Ti)     | [7440-32-6] | 1000                | mg/L |        | (NH <sub>4</sub> ) <sub>2</sub> TiF <sub>6</sub><br>,, 99.998% |