

Yttrium (Y) 1000 mg/L in 2-5% HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003017
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Yttrium (Y)	[7440-65-5]	1000	mg/L		Y2O ₃ , 99.9 99%