

## **Arsenic (As) 100 mg/L in HNO<sub>3</sub>, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4003021  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

### Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100	mg/L		As, 99.999 %