

Barium (Ba) 100 mg/L in diluted HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003024
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	100	mg/L		BaCO ₃ , 99.999%