

**Beryllium Be - 100 mg/L in diluted HNO3 for ICP-MS (Be4O(C2H3O2)6, 99.999%)**

Art. ID LM24-OP-4003026  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Beryllium (Be)    | [7440-41-7] | 100                 | mg/L |        |        |