

Beryllium Be - 100 mg/L in diluted HNO3 for ICP-MS (Be4O(C2H3O2)6, 99.999%)

Art. ID LM24-OP-4003026
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	100	mg/L		