

Dysprosium (Dy) 100 mg/L in diluted HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003035
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dysprosium (Dy)	[7429-91-6]	100	mg/L		Dy ₂ O ₃ , 99.999%