

Iron (Fe) 100 mg/L in 2% HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003038
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iron (Fe)	[7439-89-6]	100	mg/L		Fe(NO ₃) ₃ , 99.9995%