

**Germanium (Ge) 100 mg/L in diluted HNO3/HF tr, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4003041  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Germanium (Ge)	[7440-56-4]	100	mg/L		Ge, 99.999 %