

Hafnium (Hf) 100 mg/L in diluted HNO₃/HF tr, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003042
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Hafnium (Hf)	[7440-58-6]	100	mg/L		HfO ₂ , 99.9 95%