

Iridium (Ir) 100 mg/L in diluted HCl, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003046
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iridium (Ir)	[7439-88-5]	100	mg/L		(NH ₄) ₂ IrCl 6, 99.99%