

## **Manganese (Mn) 100 mg/L in 2-5% HNO<sub>3</sub>, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4003052  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

### Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Manganese (Mn)	[7439-96-5]	100	mg/L		Mn(NO <sub>3</sub> ) <sub>2</sub> , 99.999%