

Nickel (Ni) 100 mg/L in 2-5% HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003057
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|------------------|
| | Nickel (Ni) | [7440-02-0] | 100 | mg/L | | Ni, 99.999 8% |