

Lead (Pb) 100 mg/L in 2 % HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003059
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lead (Pb)	[7439-92-1]	100	mg/L		Pb(NO ₃) ₂ , 99.9995%