

**Titanium (Ti) 100 mg/L in diluted HNO<sub>3</sub>/HF tr, Single-Element Standard for ICP/MS**

Art. ID LM24-OP-4003078  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Titanium (Ti)	[7440-32-6]	100	mg/L		(NH <sub>4</sub> ) <sub>2</sub> TiF <sub>6</sub> ,, 99.998%