

Beryllium (Be) 10 mg/L in 2% HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003096
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	10	mg/L		Be4O(C ₂ H ₃ O 2) ₆ , 99.99 9%