

Sodium (Na) 10 mg/L in 2-5% HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003125
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|-----------------------------|
| | Sodium (Na) | [7440-23-5] | 10 | mg/L | | NaNO ₃ , 99.999% |