

Tin (Sn) 10 mg/L in diluted HNO₃/HF tr, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003144
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	10	mg/L		Sn, 99.999 9%