

Titanium (Ti) 10 mg/L in 2-5% HNO₃/HF tr, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003149
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--|
| | Titanium (Ti) | [7440-32-6] | 10 | mg/L | | (NH ₄) ₂ TiF ₆ ,, 99.998% |