

Zinc (Zn) 10 mg/L in 2% HNO₃, Single-Element Standard for ICP/MS

Art. ID LM24-OP-4003157
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Zinc (Zn)	[7440-66-6]	10	mg/L		Zn, 99.999 %