

**Phosphate Modifier for Graphite Furnace - NH<sub>4</sub>H<sub>2</sub>PO<sub>4</sub> 20 g/L in 1 %HNO<sub>3</sub>**

Art. ID LM24-OP-4003190  
Unit 50 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	NH <sub>4</sub> H <sub>2</sub> PO <sub>4</sub>		20	g/L		