

## **ICP-MS MIX 22 components in 2% HNO3**

Art. ID LM24-OP-4003192  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

### Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	5	mg/L		
	Aluminium (Al)	[7429-90-5]	30	mg/L		
	Arsenic (As)	[7440-38-2]	10	mg/L		
	Boron (B)	[7440-42-8]	30	mg/L		
	Barium (Ba)	[7440-39-3]	10	mg/L		
	Beryllium (Be)	[7440-41-7]	3	mg/L		
	Cadmium (Cd)	[7440-43-9]	3	mg/L		
	Cobalt (Co)	[7440-48-4]	5	mg/L		
	Chromium (Cr)	[7440-47-3]	10	mg/L		
	Copper (Cu)	[7440-50-8]	10	mg/L		
	Manganese (Mn)	[7439-96-5]	5	mg/L		
	Molybdenum (Mo)	[7439-98-7]	3	mg/L		
	Nickel (Ni)	[7440-02-0]	10	mg/L		
	Lead (Pb)	[7439-92-1]	5	mg/L		
	Selenium (Se)	[7782-49-2]	5	mg/L		
	Tin (Sn)	[7440-31-5]	10	mg/L		
	Thallium (Tl)	[7440-28-0]	2	mg/L		
	Tellurium (Te)	[13494-80-9]	10	mg/L		
	Vanadium (V)	[7440-62-2]	10	mg/L		
	Zinc (Zn)	[7440-66-6]	100	mg/L		
	Strontium (Sr)	[7440-24-6]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	3	mg/L		