

ICP/MS Calibration standard 5 components in Nitric Acid 2%

Art. ID LM24-OP-4003196
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	1000	mg/L		
	Magnesium (Mg)	[7439-95-4]	1000	mg/L		
	Potassium (K)	[7440-09-7]	1000	mg/L		
	Sodium (Na)	[7440-23-5]	1000	mg/L		
	Iron (Fe)	[7439-89-6]	1000	mg/L		