

ICP Multi-Element Standard - Tuning Solution 4 elements in 2 %HNO3

Art. ID LM24-OP-4003201
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cerium (Ce)	[7440-45-1]	10	mg/L		
	Lithium (Li)	[7439-93-2]	10	mg/L		
	Thallium (Tl)	[7440-28-0]	10	mg/L		
	Yttrium (Y)	[7440-65-5]	10	mg/L		