

ICP Multi-Element Standard - Quality Control-MS Multi 8 elements in 2 %HNO3/HF tr

Art. ID LM24-OP-4003225
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Germanium (Ge)	[7440-56-4]	10	mg/L		
	Hafnium (Hf)	[7440-58-6]	10	mg/L		
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Tellurium (Te)	[13494-80-9]	10	mg/L		
	Tin (Sn)	[7440-31-5]	10	mg/L		
	Tungsten (W)	[7440-33-7]	10	mg/L		
	Zirconium (Zr)	[7440-67-7]	10	mg/L		