

ICP-MS MIX 16 components in 5% HNO3

Art. ID LM24-OP-4003242
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	10	ug/ml		
	Strontium (Sr)	[7440-24-6]	10	ug/ml		
	Manganese (Mn)	[7439-96-5]	10	ug/ml		
	Lithium (Li)	[7439-93-2]	10	ug/ml		
	Cobalt (Co)	[7440-48-4]	10	ug/ml		
	Arsenic (As)	[7440-38-2]	10	ug/ml		
	Zinc (Zn)	[7440-66-6]	10	ug/ml		
	Selenium (Se)	[7782-49-2]	10	ug/ml		
	Caesium (Cs)	[7440-46-2]	10	ug/ml		
	Rubidium (Rb)	[7440-17-7]	10	ug/ml		
	Lead (Pb)	[7439-92-1]	10	ug/ml		
	Scandium (Sc)	[7440-20-2]	10	ug/ml		
	Nickel (Ni)	[7440-02-0]	10	ug/ml		
	Vanadium (V)	[7440-62-2]	10	ug/ml		
	Chromium (Cr)	[7440-47-3]	10	ug/ml		
	Beryllium (Be)	[7440-41-7]	10	ug/ml		