

## **UV Wavecal Solution - 12 components (Standard Equivalent to Perkin Elmer)**

Art. ID LM24-OP-4003249

Unit 500 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphorus (P)	[7723-14-0]	100	mg/L		
	Sulfur (S)	[7704-34-9]	100	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	20	mg/L		
	Molybdenum (Mo)	[7439-98-7]	20	mg/L		
	Nickel (Ni)	[7440-02-0]	20	mg/L		
	Scandium (Sc)	[7440-20-2]	20	mg/L		
	Arsenic (As)	[7440-38-2]	20	mg/L		
	Sodium (Na)	[7440-23-5]	20	mg/L		
	Lanthanum (La)	[7439-91-0]	20	mg/L		
	Lithium (Li)	[7439-93-2]	20	mg/L		
	Calcium (Ca)	[7440-70-2]	1	mg/L		