

Mixed Calibration Standard - 10 components (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003253
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	50	mg/L		
	Potassium (K)	[7440-09-7]	50	mg/L		
	Lanthanum (La)	[7439-91-0]	10	mg/L		
	Lithium (Li)	[7439-93-2]	10	mg/L		
	Manganese (Mn)	[7439-96-5]	10	mg/L		
	Nickel (Ni)	[7440-02-0]	10	mg/L		
	Strontium (Sr)	[7440-24-6]	10	mg/L		
	Zinc (Zn)	[7440-66-6]	10	mg/L		
	Barium (Ba)	[7440-39-3]	1	mg/L		
	Magnesium (Mg)	[7439-95-4]	1	mg/L		