

Low UV Standard - 3 components (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003257
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	10	mg/L		
	Phosphorus (P)	[7723-14-0]	10	mg/L		
	Sulfur (S)	[7704-34-9]	10	mg/L		