

Elan 5000 Detection Limit Solution - 6 components 10 µg/L each in HNO3 2% (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003261
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 12 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Beryllium (Be) | [7440-41-7] | 10 | µg/L | | |
| | Cobalt (Co) | [7440-48-4] | 10 | µg/L | | |
| | Germanium (Ge) | [7440-56-4] | 10 | µg/L | | |
| | Indium (In) | [7440-74-6] | 10 | µg/L | | |
| | Thallium (Tl) | [7440-28-0] | 10 | µg/L | | |
| | Uranium (U) | [7440-61-1] | 10 | µg/L | | |