

Elan 5000 Detection Limit Solution - 6 components 10 µg/L each in HNO₃ 2% (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003262
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	10	µg/L		
	Cobalt (Co)	[7440-48-4]	10	µg/L		
	Germanium (Ge)	[7440-56-4]	10	µg/L		
	Indium (In)	[7440-74-6]	10	µg/L		
	Thallium (Tl)	[7440-28-0]	10	µg/L		
	Uranium (U)	[7440-61-1]	10	µg/L		