

Elan 9000/6100 Setup/Stability/Masscal Solution -9 components (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003264
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 12 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Magnesium (Mg) | [7439-95-4] | 10 | µg/L | | |
| | Copper (Cu) | [7440-50-8] | 10 | µg/L | | |
| | Rhodium (Rh) | [7440-16-6] | 10 | µg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 10 | µg/L | | |
| | Indium (In) | [7440-74-6] | 10 | µg/L | | |
| | Barium (Ba) | [7440-39-3] | 10 | µg/L | | |
| | Cerium (Ce) | [7440-45-1] | 10 | µg/L | | |
| | Lead (Pb) | [7439-92-1] | 10 | µg/L | | |
| | Uranium (U) | [7440-61-1] | 10 | µg/L | | |