

Elan 6100 DRC Setup/Stab/Masscal Solution -12 components (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003269
Unit 100 mL
Deliverydetails No Dangerous Good

Description

Shelf-life: 12 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Barium (Ba) | [7440-39-3] | 10 | µg/L | | |
| | Aluminium (Al) | [7429-90-5] | 1 | µg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 1 | µg/L | | |
| | Cerium (Ce) | [7440-45-1] | 1 | µg/L | | |
| | Chromium (Cr) | [7440-47-3] | 1 | µg/L | | |
| | Copper (Cu) | [7440-50-8] | 1 | µg/L | | |
| | Indium (In) | [7440-74-6] | 1 | µg/L | | |
| | Lead (Pb) | [7439-92-1] | 1 | µg/L | | |
| | Magnesium (Mg) | [7439-95-4] | 1 | µg/L | | |
| | Manganese (Mn) | [7439-96-5] | 1 | µg/L | | |
| | Rhodium (Rh) | [7440-16-6] | 1 | µg/L | | |
| | Thorium (Th) | [7440-29-1] | 1 | µg/L | | |