

Elan 6100 DRC Setup/Stab/Masscal Solution -12 components (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003270
Unit 500 mL
Deliverydetails No Dangerous Good /not restricted

Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	10	µg/L		
	Aluminium (Al)	[7429-90-5]	1	µg/L		
	Cadmium (Cd)	[7440-43-9]	1	µg/L		
	Cerium (Ce)	[7440-45-1]	1	µg/L		
	Chromium (Cr)	[7440-47-3]	1	µg/L		
	Copper (Cu)	[7440-50-8]	1	µg/L		
	Indium (In)	[7440-74-6]	1	µg/L		
	Lead (Pb)	[7439-92-1]	1	µg/L		
	Magnesium (Mg)	[7439-95-4]	1	µg/L		
	Manganese (Mn)	[7439-96-5]	1	µg/L		
	Rhodium (Rh)	[7440-16-6]	1	µg/L		
	Thorium (Th)	[7440-29-1]	1	µg/L		