

SmartTune Solution for Standard ELANs/DRC-e - 9 components (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003273
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	10	µg/L		
	Magnesium (Mg)	[7439-95-4]	10	µg/L		
	Cobalt (Co)	[7440-48-4]	10	µg/L		
	Indium (In)	[7440-74-6]	10	µg/L		
	Rhodium (Rh)	[7440-16-6]	10	µg/L		
	Cerium (Ce)	[7440-45-1]	10	µg/L		
	Barium (Ba)	[7440-39-3]	10	µg/L		
	Lead (Pb)	[7439-92-1]	10	µg/L		
	Uranium (U)	[7440-61-1]	10	µg/L		