

SmartTune Solution for DRC/DRCplus/DRC II - 10 components in HNO3 0.5% (Standard equivalent to Perkin Elmer)

Art. ID LM24-OP-4003275
Unit 500 mL
Delivery details ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) /not restricted

Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	10	µg/L		
	Beryllium (Be)	[7440-41-7]	1	µg/L		
	Cerium (Ce)	[7440-45-1]	1	µg/L		
	Cobalt (Co)	[7440-48-4]	1	µg/L		
	Indium (In)	[7440-74-6]	1	µg/L		
	Iron (Fe)	[7439-89-6]	1	µg/L		
	Lead (Pb)	[7439-92-1]	1	µg/L		
	Magnesium (Mg)	[7439-95-4]	1	µg/L		
	Thorium (Th)	[7440-29-1]	1	µg/L		
	Uranium (U)	[7440-61-1]	1	µg/L		