

Instrument Calibration Standard 1 - 4 components in HNO3 5% (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003306
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	5000	mg/L		
	Potassium (K)	[7440-09-7]	5000	mg/L		
	Sodium (Na)	[7440-23-5]	5000	mg/L		
	Magnesium (Mg)	[7439-95-4]	5000	mg/L		