

TCLP Standard 1 - 8 components (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003332
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	500	mg/L		
	Silver (Ag)	[7440-22-4]	25	mg/L		
	Arsenic (As)	[7440-38-2]	25	mg/L		
	Chromium (Cr)	[7440-47-3]	25	mg/L		
	Lead (Pb)	[7439-92-1]	25	mg/L		
	Mercury (Hg)	[7439-97-6]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	5	mg/L		
	Selenium (Se)	[7782-49-2]	5	mg/L		