

## **PerkinElmer Pure I - 19 components (Standard Equivalent to Perkin Elmer)**

Art. ID LM24-OP-4003370  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

### Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Thallium (Tl)	[7440-28-0]	400	mg/L		
	Bismuth (Bi)	[7440-69-9]	200	mg/L		
	Indium (In)	[7440-74-6]	200	mg/L		
	Lead (Pb)	[7439-92-1]	200	mg/L		
	Gallium (Ga)	[7440-55-3]	150	mg/L		
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Nickel (Ni)	[7440-02-0]	50	mg/L		
	Silver (Ag)	[7440-22-4]	50	mg/L		
	Chromium (Cr)	[7440-47-3]	25	mg/L		
	Cadmium (Cd)	[7440-43-9]	20	mg/L		
	Cobalt (Co)	[7440-48-4]	20	mg/L		
	Copper (Cu)	[7440-50-8]	20	mg/L		
	Zinc (Zn)	[7440-66-6]	20	mg/L		
	Boron (B)	[7440-42-8]	15	mg/L		
	Iron (Fe)	[7439-89-6]	15	mg/L		
	Barium (Ba)	[7440-39-3]	5	mg/L		
	Manganese (Mn)	[7439-96-5]	5	mg/L		
	Beryllium (Be)	[7440-41-7]	1	mg/L		
	Strontium (Sr)	[7440-24-6]	1	mg/L		