

*****please see LM24-CUS-73001 *** Multi-Element Solution - 5 components (Standard Equivalent to Perkin Elmer)**

Art. ID LM24-OP-4003389
Unit 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		
	Antimony (Sb)	[7440-36-0]	1000	mg/L		
	Tungsten (W)	[7440-33-7]	1000	mg/L		
	Zirconium (Zr)	[7440-67-7]	1000	mg/L		