

Multi-Element Solution - 5 components (Standard Equivalent to Perkin Elmer)

Art. ID LM24-OP-4003390
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Molybdenum (Mo) | [7439-98-7] | 1000 | mg/L | | |
| | Tin (Sn) | [7440-31-5] | 1000 | mg/L | | |
| | Antimony (Sb) | [7440-36-0] | 1000 | mg/L | | |
| | Tungsten (W) | [7440-33-7] | 1000 | mg/L | | |
| | Zirconium (Zr) | [7440-67-7] | 1000 | mg/L | | |