

## **Multi-Element Solution - 17 components (Standard Equivalent to Perkin Elmer)**

Art. ID LM24-OP-4003391  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	1000	mg/L		
	Barium (Ba)	[7440-39-3]	1000	mg/L		
	Beryllium (Be)	[7440-41-7]	1000	mg/L		
	Chromium (Cr)	[7440-47-3]	1000	mg/L		
	Cadmium (Cd)	[7440-43-9]	1000	mg/L		
	Cobalt (Co)	[7440-48-4]	1000	mg/L		
	Copper (Cu)	[7440-50-8]	1000	mg/L		
	Lanthanum (La)	[7439-91-0]	1000	mg/L		
	Lead (Pb)	[7439-92-1]	1000	mg/L		
	Lithium (Li)	[7439-93-2]	1000	mg/L		
	Manganese (Mn)	[7439-96-5]	1000	mg/L		
	Nickel (Ni)	[7440-02-0]	1000	mg/L		
	Scandium (Sc)	[7440-20-2]	1000	mg/L		
	Strontium (Sr)	[7440-24-6]	1000	mg/L		
	Vanadium (V)	[7440-62-2]	1000	mg/L		
	Yttrium (Y)	[7440-65-5]	1000	mg/L		
	Zinc (Zn)	[7440-66-6]	1000	mg/L		