

## **2,4,6-Trichloroaniline 100 µg/mL in Acetone**

Art. ID LM24-OP-4003414  
Unit 1 mL  
Deliverydetails De Minimis UN1090 II E2

### Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter      | CAS number | Concentration/Value | Unit  | Method | Source |
|------------------|------------------------|------------|---------------------|-------|--------|--------|
|                  | 2,4,6-Trichloroaniline | [634-93-5] | 100                 | µg/mL |        |        |