

## **Piperonyl butoxide 100 µg/mL in Methanol**

Art. ID LM24-OP-4005745  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

### Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter  | CAS number | Concentration/Value | Unit  | Method | Source |
|------------------|--------------------|------------|---------------------|-------|--------|--------|
|                  | Piperonyl butoxide | [51-03-6]  | 100                 | µg/mL |        |        |