

Stainless steel Z5CNU15.05 with 14,7 % Cr, 4,8 % Ni and 3 % - chips

Art. ID ECRM-F 273-1
Unit 100 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,0336	%		
	Silicon (Si)	[7440-21-3]	0,378	%		
	Manganese (Mn)	[7439-96-5]	0,785	%		
	Phosphorus (P)	[7723-14-0]	0,0131	%		
	Sulfur (S)	[7704-34-9]	0,00037	%		
	Chromium (Cr)	[7440-47-3]	14,747	%		
	Molybdenum (Mo)	[7439-98-7]	0,246	%		
	Nickel (Ni)	[7440-02-0]	4,852	%		
	Arsenic (As)	[7440-38-2]	0,003	%		
	Cobalt (Co)	[7440-48-4]	0,0391	%		
	Copper (Cu)	[7440-50-8]	3,047	%		
	Nitrogen (N)	[7727-37-9]	0,0444	%		
	Niobium (Nb)	[7440-03-1]	0,221	%		
	Tin (Sn)	[7440-31-5]	0,0021	%		
	Vanadium (V)	[7440-62-2]	0,000512	%		