

**p-Terphenyl D14 2000 µg/mL in Dichloromethane**

Art. ID LM24-OP-4007051  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1593 III E1

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
|                  | p-Terphenyl D14   | [1718-51-0] | 2000                | µg/mL |        |        |