

### **4-Chlorobiphenyl (PCB 3) 100 µg/ml in Iso-octane**

Art. ID                      LM24-OP-4007063  
Unit                         1 mL  
Deliverydetails            De Minimis UN1262 II E2

#### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Chlorobiphenyl (PCB 3 )	[2051-62-9]	100	µg/ml		